

Collected From: _____ On: _____

Preparation Method / Extender Used:

Shipping Method:

Cooled at .7 °C to 18 °C, then cooled at .05 °C per minute to 5 °C

Ship-Mate Electronic Shipping Thermos information can be found at www.itlpros.com/shipmate

The semen is inside the Ship-Mate Thermos

(We recommend AI-ing as soon as possible)

The green light on thermos lid will be flashing the temperature

5 flashes is 5 °C, which is the desired temperature.

(anywhere from 2-10 flashes is in range)

Please Return the Ship-Mate Thermos

Use the Prepaid Shipping Label

VERY IMPORTANT: Discard the Battery Before Shipping Back

(We've had a bomb squad evacuated a building because the batteries were left in and some on thought it was suspicious)

Extra Information:
